

### P-025 - OPTIMIZATION OF LINOLEIC ACID EMULSION PREPARATION TO REDUCE SUBSTRATE LOSSES AFTER FILTER-STERILIZATION

Ana Luiza Fontes<sup>1</sup>; Lúgia Pimentel<sup>1,2,3</sup>; Ana Sofia Salsinha<sup>1</sup>; Beatriz Cardoso<sup>1,4</sup>; José Carlos Andrade<sup>4</sup>; Ana Maria Gomes<sup>1</sup>; Luis Miguel Rodríguez-Alcalá<sup>1</sup>

1 - Universidade Católica Portuguesa, CBQF - Centro de Biotecnologia e Química Fina – Laboratório Associado, Escola Superior de Biotecnologia; 2 - CINTESIS – Centro de Investigação em Tecnologias e Sistemas de Informação em Saúde, Faculdade de Medicina da Universidade do Porto; 3 - QOPNA – Unidade de Investigação de Química Orgânica, Produtos Naturais e Agroalimentares, Universidade de Aveiro; 4 - Centro de Investigação em Ciências da Saúde (CICS), Instituto Superior de Ciências da Saúde – Norte, CESPU

#### Background

Conjugated linoleic acid (CLA) isomers are bioactive fatty acids that can be produced microbiologically from linoleic acid (LA). *In vitro* studies normally test CLA production using a solution of pure LA at a specific concentration and an emulsifier, usually Tween 80, upon which the mixture is filter-sterilized [1]. However, preparation leads to LA losses requiring higher amounts of LA to achieve the intended concentration. Thus, the aim of this work was to optimize the LA emulsification strategy in order to obtain a more efficient and cost-effective procedure.

#### Method

Four different treatments were applied to LA solutions at 15 mg/mL with 2% (w/v) Tween 80: i) Filtration (0.45 µm-pore size membrane; F4); ii) Ultra-Turrax (90s [A] or 150s [B] at intervals of 30s) + Filtration (0.45 µm); iii) Sonicator (90s at intervals of 30s; 1s pulses [C]) + Filtration (0.45 µm). The application of a smaller pore sized membrane assuring microbiological sterility was further tested and two different treatments were applied: i) Filtration (0.20 µm-pore size membrane; F2); ii) Ultra-Turrax (150s at intervals of 30s [D]) + Filtration (0.20 µm). All experiments were carried out in duplicate. Aliquots of each solution were collected before and after treatment for LA concentration analysis by gas chromatography [2].

#### Results & Conclusions

Filtration of LA emulsion directly through a 0.45 µm-pore size membrane (F4) led to a 17.09% LA loss whereas the introduction of a blending step, either using a sonicator (C) or an Ultra-Turrax (A and B) led to lower losses of 13.01% and 7.17% and 7.40% ( $p > 0.05$ ), respectively. A smaller filter pore size membrane (0.20 µm), also contributed to lower losses; LA reductions were from 10.07% (F2) to 3.71% (D). In conclusion, a previous dispersion with Ultra-Turrax, independent of the filter pore size, demonstrated to be the best method to reduce substrate losses in filter-sterilization of LA emulsions.

#### References & Acknowledgments

[1] A.L. Fontes, L.L. Pimentel, C.D. Simões, A.M.P. Gomes, L.M. Rodríguez-Alcalá, Crit. Rev. Food Sci. Nutr. 57 (2017) 2611–2622.

[2] L.L. Pimentel, A.L. Fontes, A.M. Gomes, L.M. Rodríguez-Alcalá, MethodsX 2 (2015) 475–484.

This work was financed by National funds via FCT – Fundação para a Ciência e a Tecnologia, under the project “Pro-TECh-CLnA - Microbial Production of Bioactive Conjugated Linolenic Acid Isomers to Obtain Functional Ingredients and Foods” reference PTDC/AGR-TEC/2125/2014. Financial support for the authors A.L. Fontes and L.L. Pimentel was provided by fellowships SFRH/BD/117721/2016 and SFRH/BPD/119785/2016, respectively, granted by the Portuguese government through FCT.

**Keywords:** Linoleic acid, Emulsions, Sterilization, CLA production substrate, Homogenization