

# Smart sampling procedure for metal ions assessment in dynamic water systems



CATOLICA

CBQF · CENTRO DE BIOTECNOLOGIA  
E QUÍMICA FINA LABORATÓRIO ASSOCIADO

CBQF

Raquel B. R. Mesquita, Mafalda G. Pereira, Francisco A. P. Rodrigues, Maria M. P. Melo, Francisca T. S. M. Ferreira, Tânia C. F. Ribas, António O. S. S. Rangel

Universidade Católica Portuguesa, CBQF - Centro de Biotecnologia e Química Fina – Laboratório Associado, Escola Superior de Biotecnologia, Porto, Portugal

rmesquita@ucp.pt

PORTO

XXVIII ENCONTRO

GALEGO  
PORTUGUÊS  
DE QUÍMICA

## Introduction

### Water bodies are dynamic systems



High necessity to monitor metal ions

Can cause adverse effects to living microorganisms:



Toxic



Carcinogenic



Mutagenic



Teratogenic

## Aqua\_Smart



### Project

Smart sampling    Real-time monitoring    In-situ systems



## GREEN CHEMISTRY



Promoting **green chemistry** by employing techniques that significantly reduce the use of hazardous substances.

## ONE HEALTH



Monitoring environmental changes in waters namely metal ions content, enables to prevent threats to public health.




## Sampling *in situ* measurements

*In situ* probe measurements:


- pH
- Conductivity
  - Salinity
- Temperature
- Suspended solids
- Dissolved oxygen

## Cadmium determination

### *In situ* Smart Sampling



### Elution and Analysis



### Main Results

- Sampling and pre-treatment performed *in situ*;
- Matrix discarded;
- Achieved a five-time pre-concentration of the sample.

## Iron(II) & iron(III) speciation

*In situ*:

Fe<sup>2+</sup> from the sample (direct determination)  
 Fe<sup>3+</sup> from the sample (bead injection)

	Dynamic Range (mg/L)	Calibration Curve
Direct determination	0.05 – 2.0	A = 0.0417 × [Fe] + 0.013
Bead injection	0.25 – 1.0	A = 0.0237 × [Fe <sup>3+</sup> ] + 0.017

**On going...**

## Cu<sup>2+</sup>, Mn<sup>2+</sup> and Zn<sup>2+</sup> determination

### Strategy




### So far...

- ✓ Retention of Cu<sup>2+</sup> using anionic exchange resin;
- ✓ Metal ions quantification with PAR;

### On going...

- ✓ Optimization of quantification with Zincon.

### System Manifold



## Acknowledgments

This work was supported by National Funds from FCT - Fundação para a Ciência e a Tecnologia through project 2022.08713.PTDC; we would also like to thank the scientific collaboration under the FCT project UIDB/50016/2020.

